# Tiger Dental

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
## Notizen

**English**

**Dear doctor,** we are honored that you have chosen this innovative new product from Tiger Dental GmbH.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at [tigerdental.com](http://tigerdental.com) under the “Know-how” link. Our customer service is available for any additional questions.

The continual development of our products indicates our company’s commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the “Downloads” link at [tigerdental.com](http://tigerdental.com)).

## **1. Manufacturer**

Tiger Dental GmbH, Belruptstraße 59, 6900 Bregenz, Austria

## **2. Product Description**

The Flex Developer (FD) is a highly efficient class-II-corrector that achieves tooth movement as well as triggers mandibular growth. Therapeutic force values remain constant even after millions of bending cycles, allowing for class-II-correction in the shortest possible time.

An additional indication is mandibular space-closure, especially following extraction or due to congenitally missing teeth.

If desired, distalization of maxillary teeth is also possible. Operators love the Flex Developer due to its robustness and fracture resistance.

In addition to moving teeth dentally, the device also has a functional effect. This stimulates the patient to position the mandible forward and induces Herbst-like mandibular growth. Biting in centric occlusion remains possible due to the device’s exceptional flexibility.

The spring action is created using a fracture resistant, constantly flexible polyamide-rod. As a result, the FD can be shortened to any length and thinned down to any diameter. In this fashion forces between 50 – 1,000 cN (cN ~ g) can be selected.

Because the FD is secured to a sectional mandibular guide wire, it delivers a mesializing force to the lower molars that is especially useful in

closing mandibular spaces.

The FD leads to bodily tooth movement of 1 mm/month and can also be used unilaterally.

### **3. Applications**

#### **Indications:**

- Distalization of maxillary molars (headgear effect)
- Mesialization of mandibular molars (during extraction treatment or congenitally missing lower premolars)
- Mandibular growth acceleration (activator effect)
- Whole arch class-II-correction (FS/removable)
- Retrusion of the maxillary anteriors (used in conjunction with elastic chain)
- Midline correction (one-sided use)
- Unilateral dental class-II-correction (one-sided use)

#### **Contraindications:**

- Protrusive mandibular anterior teeth (except when closing space in buccal segments, congenitally missing premolars)
- Steep occlusal plane (gummy smile), exercise caution with anterior open bites

## 4. Installation and removal of the Flex Developer

### System requirements

- Molar bands with traditional and headgear tubes, mandibular bands with double tubes (.018" or .022")
- Full-fixed maxillary and mandibular brackets

### Armamentarium

- FD starter kit:
  - FD 2.0 box: 5 pairs of variable length spring rods
  - FD 2.0 Pin pliers for pressing in the fixation point or clamp pliers for clamping the closing piece
  - Torque key
  - FD 2.0 Measuring ruler
- Three-prong plier (extra fine)
- Weingart plier
- Distal end cutter

### Preparing the guide wire

Fig. 2: Before installation please bend the lower molar hook anteriorly and gingivally. This way the spring rod will not become locked underneath the molar attachment causing it to loosen.

**Notice:** Right-side guide wires are indicated with red elastics and left-side guide wires are indicated with yellow elastics. Please remove these elastics before installation.


Fig. 3: Cut the guide wire

Shorten the guide wire so that it fits in the patient's mouth. Insert the guide wire in the auxiliary tube as shown until the guide wire hook is approximately half-way between the first lower premolar and the canine brackets. Curve the guide wire slightly toward the tongue so that its hook sits passively underneath the base archwire.

Fig. 4: Shorten with distal end cutter, step down bend.

Subsequently shorten the guide wire so that approximately 5 mm of the wire remain distal to the auxiliary tube. Make two marks on the wire between the first molar and second premolar as well as below the second premolar bracket.

Fig. 5 & 6: Make the first bend with the FD-three prong plier at the posterior mark. The single prong should be facing gingivally. Execute a 40 bend.

2


3


4


5

6

7

8

9


10


11


12


13

14


15

16

17a

b

Fig. 7 & 8: Make a second bend with the three prong plier flipped at the second mark. Bend until the guide wire is again parallel with the base archwire. It is important that there is sufficient separation between the guide wire step bend and the lower second premolar bracket. The FD can become wedged between the bracket and guide wire. Ultimately, the guide wire should run approximately 3-4 mm apical to the brackets.

Fig. 9 & 10: If necessary, remove the guide wire and make adjustments with your fingers to ensure a the guide wire hook's passive placement between the lower canine and first premolar (Fig. 9a: no bend placed, Fig. 9b: the curvature of the base archwire is matched).

**Notice:** Never use a rough plier to adjust the glide plane of the guide wire.

Fig. 11 & 12: Secure the guide wire hook to the base archwire by pressing it closed from the occlusal aspect with a fine Wiengart plier. This forms a ring that glides loosely on the base archwire. Tip: Use light-curing resin on the locked hook to minimize vertical play at this junction. Gliding action will be maintained.

### **The next steps depend on intended use:**

Fig. 13: Pull the guide wire distally through the auxiliary molar tube and **bend occlusally**. This way a mandibular molar or the entire buccal segment can be mesialized to close space by protracting posterior teeth.

Fig. 14: **Do not cinch the guide wire distal** to the molar tube to allow the guide wire hook to rest against the lower canine bracket. This will increase anterior tooth protrusion and will open space. In this case, **also do not cinch the lower base archwire**.

After making the decision whether or not to cinch the lower archwires, you must also decide whether to cinch the maxillary archwire. When the maxillary archwire is not cinched, the force from the spring rod will distalize the upper molar in a headgear-like fashion. This creates maxillary arch length. If the maxillary archwire is cinched, the maxilla will be distalized en-masse or can help serve as anchorage.

### **Installation of the Flex Developer:**

Fig. 15: To install the FD, you must measure from the mesial of the upper molar tube to the most anterior portion of the guide wire. Please use the side of the ruler labelled with a smiley-face 😊 (22 mm length pictured).

**Notice:** Right-side FDs are marked in red. The locking end with hook faces buccally.

### **Measuring and Aligning the Flex Developer:**

Fig. 16–19: Push the FD into the ruler's tube and align the measured distance (22 mm in this example) at the anterior portion of the looped end. Then cut off the projecting part with the sharp FD cutter. Attach the closing piece and turn towards the gingiva until the laser-marked long axis lines are aligned (Fig. 18).

Using the pin on the FD pin pliers, firmly clamp the closing piece at the fixation point (crater-shaped weak point). Alternatively clamp twice with the clamp pliers so that 2 notches are visible on the sides. (Fig. 19a+b). The hook is bent back approximately 45° from the spring rod's long axis. Do not change this inclination (Fig. 17b)!

### **Insertion of the Flex Developer**


Fig. 20–22: Using the securing disc to hold the pin (no pliers needed!), insert the pin distally through the headgear tube. Use a plier to bend the mesial portion of the pin (this secures the pin, but should not activate it). Bend the ball end (distal end) 30° apically.

Fig. 23: The patient should now move the lower jaw into a protrusive position. This way the FD can be hooked into the guide wire without tension. The Flex Developer will be activated when the patient attempts to place his/her mandible into a retrusive position. Warning: The hook should not be opened/closed more than 8–10 times (this increases fracture risk)!


### **Activation of the Flex Developer**

Fig 24 & 25: Remove the securing disc (the right-side disc is red, the left-side disc is white). Pull the mesial end of the ball pin through the headgear tube and bend up so that the ball is moved closer to the headgear tube. This will bend the spring rod (Fig. 31). The force of the locking end hook can be up to 400 cN (this can be confirmed with a Correx gauge, Fig. 30). Secure the mesial end of the pin with an Alastic ligature to a gingivally directed maxillary molar hook. This fixes the pin while allowing a short glide plane distal to the headgear tube. Additional activations are possible at any time through additional mesial movement of the ball pin.

18


19a Pin pliers, b Clamp pliers

20

21

22

23

24

25

26


27


28


29


30

31

32

33

## Fine adjustments of the guide wire and force measurement

Fig. 26–28: The locking end is activated with a Weingart plier. If necessary, the FD can be removed from the mouth every 30 days to perform routine checks and can be reinserted.

Fig. 29: **Tip**: If the guide wire and the locking end are too close to the gingiva (tissue blanching), use a torque key or torqueing pliers to place buccal torque in the guide wire mesial to the molar tube. Never adjust the torque in the glide plane portion of the guide wire!

Fig. 30: The force of the Flex Developer can be set individually and measured with a Correx gauge. The measuring direction is toward the pin ball. Push the gauge against the locking end hook until it begins to slide on the guide wire. The force generated when the mouth is closed is the force potential of the FD.

## Reactivation of the Flex Developer

Fig. 31–33: There are 3 possibilities to reactivate the FD and increase the force potential: **1.** Shorten the pin ball. **2.** Pull the guide wire distally. The guide wire is moved distally through the molar tube and any distal excess can once again be bent occlusally. It is possible that the step bend will need to be moved mesially. **3.** The guide wire's glide plane can be shortened using a ball of light activated composite (Heliosit®). Form a ball of composite around the guide wire with an explorer and light cure it. The guide wire eyelet must retain sliding ability!

## Removal of the Flex Developer and patient self-help

In order to remove the FD, bend open the locking hook with Weingart pliers and unhook the FD as the patient moves his/her mandible into a protrusive position. Straighten the ball hook and remove by pulling the ball distally. The patient can also temporarily deactivate the FD. To do this, a small plier is used (i.e. fingernail forceps) to slightly open the locking hook. When the patient moves his/her mandible into a protrusive position the FD will disengage and can be temporarily hooked to the maxillary base archwire.

## **5. Contents**

REF 10860 FD Flex Developer Box 2.0, variable length, 5 pairs/pack

REF 10904 FD 2.0 Flex Developer-Box, right side, 10 items/pack

REF 10905 FD 2.0 Flex Developer-Box, left side, 10 items/pack

## **6. Notices for Single-use Products**

This product is nearly ready to use out of the box, but is not sterile and should not be used more than once. The reprocessing of this once used product (recycling) and its reapplication to the patient is strictly prohibited.

## **7. Safety Instructions**

This operating and installation manual contains tips and recommendations, but does not release you from personal responsibility! The use of this product is allowed only by orthodontists and dentists according to MDD 93/42/EEC. The user must comply with the provisions of the Medical Devices Directive and any other applicable laws and regulations.

 This product contains nickel.

## **8. References of Quality**

Tiger Dental assures the user a flawless quality of its products. The content of this manual is based on experience. The user is responsible for the processing of the products themselves. Tiger Dental is not liable for erroneous results or processing by the user that occur in the absence of influence of Tiger Dental. Although this product is a mature product, improvements and changes can be made, which may not be documented in the present manual. Similarly, we are not liable for printing errors.

## 9. Description of the used symbols of ÖNORM EN ISO 15223-1

 0483 Notified body

 Manufacturer

 Date of manufacture

 Lot / Batch number

 Reference / Article number


 Non sterile

 Do not reuse

 Consult instructions for use

 Caution

## Notes

Informationen zu Produkten und Serviceleistungen finden Sie unter:  
For more information on our products and services, please visit:  
**[www.tigerdental.com](http://www.tigerdental.com)**

[www.flex-developer.com](http://www.flex-developer.com)

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