Innovations in orthodontics

Tiger Dental

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Achtung

English

Dear doctor, we are honored that you have chosen this innovative new product from Tiger Dental GmbH.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at tigerdental.com under the “Know-how” link. Our customer service is available for any additional questions.

The continual development of our products indicates our company’s commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the “Downloads” link at tigerdental.com).

1. Manufacturer

Tiger Dental GmbH, Belruptstrasse 59, 6900 Bregenz, Austria

2. Product Description

The Herbst Developer (HD) is a fixed orthodontic instrument for treatment of skeletal Class II malpositions. Mainly used between the ages of 11 and 16 years, it serves growth promotion and can be combined with a fixed functional therapy. Since the Herbst Developer corresponds to a ball joint at both ends, sideways chewing movements are also possible for the patient and the mouth opening is not hindered by the telescopic design of the instrument. The easily bent square-edged telescope can be adjusted to the desired length and its slim design and curved form make it ideal for fitting into the lateral oral vestibule.

3. Scope of application

Indications:

- Treatment of skeletal Class II malocclusions
- Promoting skeletal mandibular growth during the pre- and post-pubertal growth phase
- Closing light gaps within the mandibula
- Light Headgear effect in the maxilla
- In a limited form for therapy of mandibular joint problems

Contraindications:

- Patient too old, not recommended from the age of 25 (see literature)
- proclinated lower front teeth at the beginning of the treatment

4. Fitting of the Herbst Developer

System requirements:

- molar straps with the tube and headgear tubes in the maxilla and the mandibula with the double tubes (.18' or .22')
- fixed functional therapy in the maxilla and mandibula on a minimum of 16x22 steel arch

Basic configuration:

- Herbst Developer starter set: HD 2 sets, right sliding arch, left sliding arch, distance clips, HD forceps, HD measuring ruler
- wire pliers (super fine)
- Weingart pliers
- distal end cutter
- Orthoholder, bite forms to attach the bite blocks

Note: all right parts of the Herbst Developer (sliding arch, pin, telescopic tube, inner part) are provided with red Alastic rubbers and must be inserted on the right side of the patient. Left parts are marked with yellow Alastics. With the exception of the blue distance Alastic attached to the pin, all marking rubbers must be removed before the fixing.

Completed levelling of the dental arches with alloyed steel arches (minimum of 16x25 SS) in the upper and lower jaw. The end of the steel arches should be bent distally to the last tube in both the upper and lower jaw. Straps with headgear tubes must be placed on the first molars of the upper jaw. The straps of the first lower molars must be fitted with double tubes.

Sliding curves preparation

Fig. 2: Before fitting, the lower molar hook is to be bent to caudal and gingival. This prevents the spring stick from getting stuck under the molar attachment and possibly loosening the strap.

NB: Right sliding arches are marked with red rubbers, the left ones with yellow ones. Remove rubber before inserting.

Fig. 3: Cutting of the sliding arch In principle, the sliding arch is to be

shortened slightly so that it can be adjusted in the mouth. Insert it into the auxiliary tube until the sliding arch hook lies approx. between the first premolar and canine bracket. In addition, slightly curve the sliding arch inwards so that it runs passively below the main arch.

Fig. 4: Shorten the bypass bend with a rear distal end cutter. Then cut the sliding arch with a distal end cutter so that approx. 5 mm from the auxiliary tube protrude distally. Now make two marks on the sliding arch, one between the first mandible molars and second premolars, the other below the bracket of the second premolar.

Fig. 5, 6: Make the first bend downwards with an FD wire cutter at the rear marking (bending angle of approx. 40°).

Fig. 7, 8: With an upwards bend, bend the sliding arch parallel to the main treatment curve (bayonet or step bend) The rear premolar bracket must also show a sufficient vertical distance to the bayonet bend of the sliding arch (important, otherwise FD will get stuck here!). Finally, the sliding arch should be approx. 3 - 4 mm below the brackets.

Fig. 9, 10: If necessary, remove the sliding arch again and adjust it with the fingers, seen from above, to the dental arch so that the hook lies passively and tension-free on the main arch between the triple and the quadruple bracket (a unbent, b adapted to the course of the treatment arch).

Caution: Never use rough pliers to work on or hold the sliding arc!

Fig. 11, 12: Press the occlusal sliding arch hooks together with a pointed Weingart pliers. A ring which can slide loosely on the main arch is formed.

Tip: Light-curing resin applied to the closed hook eliminates vertical slackness. Movability is retained.

Fig. 13: Bend upwards the guide wire distally of the auxiliary molar tube. Bend the upper base archwire also upwards distally. the lower base archwire keeps straight and can be shortened distally.

2


3


4


5

6

7

8

9


10


11

12


13

14

15

16

Determination of the HD length

Fig. 14, 15: Usually the HD is correct in its length. The HD measuring ruler has two small hooks on one side, where there is an open holder with two small arms. Both sides of the ruler have a centimetre scale. One side is provided with a smiley 😊.

The patient is asked to place the 4 lower incisors behind the bite block. The measuring ruler is now attached to the back end of the headgear tube of the upper first molar with the hook. The distance from the rear edge of the headgear tube to the front kink of the sliding arch is measured.

Subsequently, the measuring ruler is placed with the smiley pointing downwards so that the two arms of the open holder point upwards. The HD is removed from the packaging and placed between the two arms of the ruler holder so that the lug and pin protrude. Then place the eyelet of the inner rod on the measuring ruler so that it points to the measured value (e.g., 27 mm).

The length control of the HD with a measuring ruler ensures that the HD protrudes approx. 5 mm distal to the headgear tube and provides place for activation. There are 2 Alastic rings in this area, which hold the lug of the telescopic tube securely pressed against the pin head. This help to avoid the annoying crackling or rubbing noises for the patient.

Fig. 16: If necessary, the length of the telescopic tube can be shortened or extended by fitting spacers. The shortening is done with a diamanted rotating instrument. The spacers are pressed onto the inner rod close to the eyelet using the HD forceps.

Inserting the HD

Fig. 17: First, straighten the hook of the molar attachment (Twister).

Fig. 18: The pin of each telescope is inserted distally into the headgear tube. Do not remove the blue distance alastics already attached to the pin. They serve as spacers.

Fig. 19: The front end of the pin is rolled up in a circle and the pin is shortened. Then the end of the pin must be fixed with the vertically bent hook of the strap attachment with an Alastic green showed in the figure. This stabilises the pin in its position.

Fig. 20: Now the inner rod is pushed in and the open eyelet at its front end can be hooked onto the sliding arch. This eyelet is now closed with a pair of Weingart pliers.

The HD of the opposite side is installed in the same way. The fine adjustment is done by shortening or extending the pin on the respective side. The aim is for the center line to coincide and the lower incisors to be directly behind the upper incisors.

Reactivation of the Herbst Developer

In order to bring the mandibula a little further forward, it may be necessary to activate the Herbst telescopes during the course of treatment. There are the following possibilities:

1. By repeated shortening the pin again (rolling of the pin end).
2. By shortening the sliding distance of the conductive arc (pulling out the bypass arc distally from the auxiliary tube and repeated bending the end of the pin).
3. Fig. 21–22: Inserting a distance clip between the kink of the inner rod and the opening of the telescopic tube. The distance clip is fixed with the HD pliers by clamping and can be removed with a bracket pliers.

Attaching palatal bite blocks

Fig. 23: To ensure secure support of the lower front teeth, bite blocks are ideally placed on the palatal side on the upper four incisors. (Mini Mold with Orthoholder).

17

18

19

20

21

22

23

5. Removal of the Herbst Developer

The telescopic tube can be removed distally by straightening the pin. The inner rod is removed by opening the eyelet with a Weingart pliers. The sliding arc can be cut and then removed.

6. Treatment period

The usual treatment period of a Herbst hinge is between 6 and 12 months depending on the severity of the malocclusion and must be determined individually for each patient depending on the course of treatment (see further literature).

7. Sources of error and difficulties

Repairs can become necessary under a high load of the Herbst hinge and unfortunately occur time and again in bruxism. The pin serves as a safety connection which can be reactivated. The pin slides out (rolls out) distally. The HD is deactivated and must be reactivated during the check-up. A fracture of the sliding curve can rarely occur under high loads. It can be readjusted and easily replaced. The eyelet of the inner rod can be opened and the new bypass arc can be hooked in and closed again.

8. Retention

After the completion of the therapy with the Herbst Developer, the bite blocks on the front teeth can remain in the mouth for several weeks. It has stood the test to have the patient wear either Class II Alastics or an activator after the removal.

9. Notices for Single-use Products

This product is nearly ready to use out of the box, but is not sterile and should not be used more than once. The reprocessing of this once used product (recycling) and its reapplication to the patient is strictly prohibited.

10. Safety instructions

This operating and installation manual contains tips and recommendations, but does not release you from personal responsibility! The use of this product is allowed only by orthodontists and dentists according to MDD 93/42/EEC. The user must comply with the provisions of the Medical Devices Directive and any other applicable laws and regulations.

 This product must be removed before magnetic resonance imaging (MRI). This product contains nickel.

11. References of Quality


Tiger Dental GmbH assures the user a flawless quality of its products. The content of this manual is based on experience. The user is responsible for the processing of the products themselves. Tiger Dental GmbH is not liable for erroneous results or processing by the user that occur in the absence of influence of Tiger Dental GmbH. Although this product is a mature product, improvements and changes can be made, which may not be documented in the present manual. Similarly, we are not liable for printing errors.

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